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Substances	Substances	Substances	Substances
11-nor- Δ -8-carboxy-THC-9-COOH (10 μ g/ml)	Deoxyepinephrine	Imipramine	Penicillin G
11-nor- Δ -9-carboxy-THC-9-COOH (10 μ g/ml)	Dextromethorphan	Isoproterenol	Phenobarbital
3,4 Methylendioxyamphetamine	d-Amphetamine	Indomethacin	Phencyclidine
4-acetaminophenol	digitoxin	Ketoprofen	Phenobarbital
Acetaminophen	digoxin	Ketamine	Phentermine
Acetylsalicylic acid	Diazepam	Lidocaine	Propanol
Amikacin	Diethylpropion	Levorphanol	Phenylethylamine-a
Amitriptyline	Diphenhydantoin	L-phenylephrine	promethazine
Ampicillin	Doxylamine	Meperidine	Quinine
Arterenol	d-Propoxyphene	Methadone	Ranitidine Sodium
Aspartame	Ecgonine HCl	Methaqualone	Salicylate
Atropine sulfate	Ecgonine Methyl Ester	Methylphenidate	Salicylic Acid
Benzoic Acid	Epinephrine(\pm)	Morphine	Secobarbital
Benzoyllecgonine	Glucose	Naloxone	Sulfate
Caffeine	Glycerol ether	Neomycin	Tetracycline
Camphor	Guaiacol	Niacinamide	Tetrahydrozoline
Chloroquine	Glucuronide	Oxazepam	Theophylline
Chlorpheniramine	Histamine	Oxycodone	Thioridazine
Chlorpromazine HCl	Homatropine	Perphenazine	Trifluoperazine
Cortisone	Hydrochlorothiazide	Phendimetrazine	Tryptophan
Cimetidine	Hydrocodone	Phencyclidine	Δ -9-THC
Codeine	Hydromorphone	Phencyclidine	

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Substances	Substances	Substances	Substances
11-nor- Δ -8-carboxy-THC-9-COOH (10 μ g/ml)	Deoxyepinephrine	Imipramine	Penicillin G
11-nor- Δ -9-carboxy-THC-9-COOH (10 μ g/ml)	Dextromethorphan	Isoproterenol	Phenobarbital
3,4 Methylendioxyamphetamine	d-Amphetamine	Indomethacin	Phencyclidine
4-acetaminophenol	digitoxin	Ketoprofen	Phenobarbital
Acetaminophen	digoxin	Ketamine	Phentermine
Acetylsalicylic acid	Diazepam	Lidocaine	Propanol
Amikacin	Diethylpropion	Levorphanol	Phenylethylamine-a
Amitriptyline	Diphenhydantoin	L-phenylephrine	promethazine
Ampicillin	Doxylamine	Meperidine	Quinine
Arterenol	d-Propoxyphene	Methadone	Ranitidine Sodium
Aspartame	Ecgonine HCl	Methaqualone	Salicylate
Atropine sulfate	Ecgonine Methyl Ester	Methylphenidate	Salicylic Acid
Benzoic Acid	Epinephrine(\pm)	Morphine	Secobarbital
Benzoylecgonine	Glucose	Naloxone	Sulfate
Caffeine	Glyceryl ether	Neomycin	Tetracycline
Camphor	Guaiacol	Niacinamide	Tetrahydrozoline
Chloroquine	Glucuronide	Oxazepam	Theophylline
Chlorpheniramine	Histamine	Oxycodone	Thioridazine
Chlorpromazine HCl	Homatropine	Perphenazine	Trifluoperazine
Cortisone	Hydrochlorothiazide	Phendimetrazine	Tryptophan
Cimetidine	Hydrocodone	Phencyclidine	Δ -9-THC
Codeine	Hydromorphone	Phencyclidine	

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- Baselt, R.C., Disposition of Toxic Drugs and Chemicals in Man, 2nd Ed., Biomedical Publ., Davis, CA, p.488(1982)
- Cody, J.T., and Schwarzhoff, R.,J. Anal. Toxicol., 17:2630 (1993)
- Urine Testing for Drugs of Abuse, NIDA Research Monograph, 73, (1986)
- Dasgupta, A., Saldana, S., Kinnaman, G., Smith, M., and Johansen, K., Clin. Chem., 39(1): 104-108 (1993)
- Department of Health and Human Services, Fed. Regist., 53(69): 11970-89 (1988)
- FDA Guidance for labeling Screening Testing, Kshitij Mohan, 7/21/1

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BIO-RAD

Liquichek™ Qualitative Urine Toxicology Control

REF	454	Negative	6 x 3 mL
	455	Positive	6 x 3 mL
	454X	MiniPak	2 x 3 mL

EXP 2020-06-30

LOT	70740	Negative	70741
		Positive	70742

ENGLISH

INTENDED USE

Liquichek Qualitative Urine Toxicology Control is intended for use as an assayed quality control urine to monitor the performance of laboratory procedures for qualitative urine toxicology.

SUMMARY AND PRINCIPLE

The use of quality control materials is indicated as an objective assessment of the precision of methods and techniques in use and is an integral part of good laboratory practices. For customers in Germany: Quality control materials are required for assessment of laboratory performance as described in the "Guideline for Quality Assurance of Medical Laboratory Examinations following the German Medical Association" (Riil-BÄK regulation).

REAGENT

This product is prepared from human urine with added drugs of abuse and metabolites of drugs of abuse, preservatives, stabilizers, and constituents of animal origin. The control is provided in liquid form for convenience.

STORAGE AND STABILITY

This product will be stable until the expiration date when stored unopened at 2 to 8°C. Once the control is opened, all analytes will be stable for 30 days when stored tightly capped at 2 to 25°C.

This product is shipped under refrigerated conditions.

PROCEDURE

This product should be treated the same as patient specimens and run in accordance with the instructions accompanying the instrument, kit, or reagent being used. Before sampling, allow the control to reach room temperature (16 to 25°C) and swirl gently to ensure homogeneity. After each use, promptly replace the stopper and return to 2 to 25°C storage.

Dispose of any discarded materials in accordance with the requirements of your local waste management authorities. In the event of damage to packaging, contact the local Bio-Rad Laboratories Sales Office or Bio-Rad Laboratories Technical Services.

LIMITATIONS

1. This product should not be used past the expiration date.
2. If there is evidence of microbial contamination or excessive turbidity in the product, discard the vial.
3. This product is not intended for use as a standard.
4. Cannabinoid metabolites may bind to certain types of plastic and glass. Refer to instructions accompanying the instrument, kit or reagent being used.
5. This product contains gentamicin. Follow the instructions provided by manufacturers of the reagent and/or test system for samples containing gentamicin.

ASSIGNMENT OF VALUES

The results printed in this insert were derived from replicate analyses and are specific for this lot of product. The tests listed were performed by the manufacturer and/or independent laboratories using manufacturer supported reagents and a representative sampling of this lot of control. Variations over time and between laboratories may be caused by differences in laboratory technique, instrumentation and reagents, or by manufacturer test method modifications. This control has been tested by gas chromatography (GC) and other methods to confirm the recovery of quantities added. These values are provided for informational use only.

Each laboratory should use the results provided only as a reference and establish its own parameters of precision.

Refer to www.qcnet.com for insert update information.

SPECIFIC PERFORMANCE CHARACTERISTICS

This product is a stabilized liquid product manufactured under rigid quality control standards. To obtain consistent vial-to-vial assay values, the control requires proper storage and handling as described.

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ENGLISH

- Qualitative results determined according to instructions from kit manufacturer.
- Some of the assays listed may not be available in all countries.
- The data required to establish the means and acceptable ranges for this assay were not obtained due to limited assignment participation. If your facility is interested in participating in the Value Assignment Program for this assay, please contact your local Bio-Rad office.

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QUALITATIVE RESULTS		
Qualitative Ergebnisse // Résultats qualitatifs // Risultati qualitativi // Resultados cualitativos // Resultados qualitativos // Kvalitativa resultat // Kvalitative resultater		
Manufacturer / Method	Negative Level ⁽¹⁾ 70741	Positive Level ⁽¹⁾ 70742
BIO-RAD TOX/SEE (2)		
Amphetamines	-	+
Barbiturates	-	+
Benzodiazepines	-	+
Cannabinoids	-	+
Cocaine	-	+
Methamphetamine	-	+
Methadone	-	+
Methylenedioxymethamphetamine (MDMA)	-	+
Opiates (300 ng/mL)	-	+
Opiates (2000 ng/mL)	-	+
Oxycodone	-	+
Phencyclidine (PCP)	-	+
Propoxyphene	-	+
Tricyclic Antidepressants (TCA)	-	+
BIOSITE TRIAGE PLUS TCA		
Amphetamines	§	§
Barbiturates	§	§
Benzodiazepines	§	§
Cannabinoids	§	§
Cocaine	§	§
Opiates	§	§
Phencyclidine (PCP)	§	§
Tricyclic Antidepressants (TCA)	§	§
BTNX RAPID RESPONSE		
d-Amphetamine	§	§
d-Methamphetamine	§	§
3, 4-Methylenedioxy-methamphetamine (MDMA)	§	§
Secobarbital	§	§
Oxazepam	§	§
11-nor-Δ-9-THC-9-COOH	§	§
Benzoylcgonine	§	§
Methadone	§	§
Methaqualone	§	§
Morphine, Free	§	§
Oxycodone	§	§
Phencyclidine (PCP)	§	§
Propoxyphene	§	§
Nortriptyline	§	§

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GPO Meth Cassette

METHAMPHETAMINE TEST (Equipment included)

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(-)-Menthhol	Berberine
(+/-)-Brom-phenil-amine	beta-Estradiol
(+/-)-Phenyl-propanol-amine	Betameth-asonone
(+/-)-Sotalol	Bilirubin
(+/-)-Suloride	Bumet-amide
(1S,9R)-beta-Hydrastine	Bupi-vacaine
(R)-(-)-Deprenyl	Butyro-phenone
1,1-Di-methyl-biguamide	Camphor
1,3-Dimethyl-uric Acid	Canrenoic Acid
1,7-Dimethyl-xanthine	Captopril
17-alpha-Ethynyl-estradiol	Carbamyl-beta-Methyl-chlorine Chloride
2-(Chloro-phenoxy)-2-Methyl-propionic Acid Ethyl Ester	Carbo-platin
2,2,2-Trichloro-ethane	Cariso-prodol
2,3-Di-mercapto-propanol (Dimer-caprol)	Cefactor
2-Butynoic Acid	Cefadroxil
2-Ethyl-2-Phenyl-malon-amide	Cefoxitin
3,3-Methyl-ene-bis-(4-Hydroxy-coumarin)(Di-cumarol)	Ceftri-axone
3,4-Dihydroxy-phenyl-acetic Acid	Cefur-oxime
3,5,5-Tri-methyl-oxazol-dine- 2,4-Dione (Trimeth-adione)	Cepha-loridine
3-Methoxy-tyramine	Cephra-dine
4,4-Diamino-diphenyl Sulfone (Dapsone)	Chlor-cyclizine
4-Dimethyl-amine-antipyrine	Chloro-thiazide
4-Hydroxy-3-Methoxy-Phenyl-acetic Acid	Chloro-trianisene
5-Hydroxy-indole-2-Carboxylic Acid	Chlor-zoxazoneZopiclone
5-Hydroxy-indole-3-Acetic Acid	Chlo-esterol
5-Hydroxy-tryptamine (Sero-tonin)	Cin-chonidine
5-Methoxy-tryptamine	Cinoxacin
5-Pregnan-3 beta-OL-20-ONE	Clobeta-sonone Butyrate
6-alpha-Methyl-17-alpha-Hydroxy- Pro-gesterone	Clozapine
6-alpha-Methyl-predni-solone	Colchicine
6-Mercap-topurine	Cortisone
Acet-aldehyde	Cromolyn
Acet-aminophen (4-acet-amido-phenol)	Cyclo-phosph-amide
Acetazo-amide	D-Amygdalin
Aceto-phenetidin	Dan-trolene
Aceto-promazine	Deferox-amine Mesylate
Albumin Standard	Desoxi-metasone
Albuterol (Sal-butamol)	Dexa-methasone
Allo-purinol	Diazoxide
alpha-Chloralose	Dichloro-methane
Alprenolol	Dichlor-phenamide
Amcino-nide	Diclofenac
Amcionide	Di-cyclomine
Amiloride	Dieldrin
Amonium Chloride	Diethyl-dithio-carbamic Acid
Ampho-tericin B	Di-florasoneDiacetate
Aniline	Diffucort-olone (pivalate)
Antipyrine	Digtoxin
Astem-izole	Dimen-hydrinate
Atenolol	Dipyrid-amide
Aza-thioprine	Dipyrene
Barbituric Acid	DL-3,4-Di-hydroxy-mandelic Acid
Beclio-methasone	DL-Aspartic Acid
Bendroflu- methiazide	DL-Indole-3-Lactic Acid
Bendro-flumethil-azide	DL-Nor-meta-nephtrine
Benzidine	DL-Phenyl-alanine
Benzic Acid beta-Diethyl-amino-ethyl Ester (Ben-actyzine)	DL-Pro-pranclol
Benzo-caine	DL-Tropic Acid
Benzoic Acid	DL-Tryptop-phanZopiclone
Benzphet-amine	d-Norprop-oxyphe
Benz-thiazide	Dobut-amine


Benzyl-amine	Doxy-cycline
d-Propoxy-phenone	Mianserin
Dro-peridol	Milrinone
d-Tubo-curarine Chloride	Minaprine
EDTA (Ethyleno-diamine-tetra-acetic Acid)	Nabu-metone
Emetine	N-Acetyl-procain-amide
Eserine	Nadolol
Estrone	Nal-buphine
Estrone-beta-D-Glucur-onide	Nalidixic Acid
Ethacrynic Acid	Nalmefene
Etham-butol	Naphthol
Ethamivan	N-Des-methyl-clozapine
Etho-propazine	Neomycin Sulfate
Ethylene Glycol	Nitro-furantoin
Ethyl para-amino-benzoate	Nor-floxacin
Etodolac	Orotic Acid
Etosopide	Oxolinic Acid
Famo-tidine	Ox-prenolol
Ferrous Sulfate	Oxy-butynin Chloride
Flu-fenamic Acid	Oxyphen-butazone
Flunisolide	Oxy-purinol
Fluphen-azine	Pacitaxel (Taxol)
Flurand-renolide	Pan-curonium Bromide
Flurbi-profen	para-Amino-benzoic Acid
Form-aldehyde	para-Phenylenediamine
Gem-fibrozil	Pargyline
Gent-amicin Sulfate	Penta-chloro-phenol
Glucose	Pent-oxifylline
Glybencl-amide (Gly-buride)	Pentyl-enetetra-zole
GuaiacolGlyceryl Ether Carbamate	Phenol
Guan-ethidine	Phenol-phthalein
Halcino-nide	Phenoxy-methyl-Penicillinic Acid (Penicillin V)
Halo-peridol	Phenol-amine
Hemo-globin	Phenyl-tolox-amine
Hexa-chloro-cyclo-hexane	Picrotoxin
Hexa-chloro-phenone	Pilo-carpine
Hippuric Acid	Pimozide
Homa-tropine	Pindolol
Hydral-azine	Pipemidic Acid
Hydrocort-tisone	Potassium Iodide
Hydro-flumethil-azide	Prazosin
Hydroxo-cobalamin	Prilocaine
Hydroxy-zine	Prima-quine
Imidazole-4-Acetic Acid	Proadifin
Indap-amide	Pro-benecid
Iptra-tropium Bromide	Pro-cyclidine
Iproniazid	Propionyl-promazine
Isoniazid Acid (Iso-nicotinic Acid Hydr-azide)	Phthalic Acid
isoprop-amide Iodide	Pyridine-2-Aldoxime- Metho-chloride
Kynurenic Acid	Pyridoxine
L-Asparic Acid	Quinolinic Acid
L-Hyoscy-amine	Rescinn-amine
Lisinopril	Reserpine
Lithium Carbonate	Ritodrine
Loper-amide	S-(-)-Carbidopa
L-Pipecolic Acid	Salicyl-amide
Mebend-azole	Salicylic Acid
Medizine	Sodium Formate
Meclo-fenamic Acid	Strychnine
Mefen-amic Acid	Succinyl- choline Chloride
Melphalan	Sulfa-methazine
Mephen-eslin	Sulfa-methi-oxazole
Mersalyl Acid	Sulfamyl-amide
Meta-camrol	Sulfathiazole
Meta-proterenol	Sulf-ioxazole
Methazol-amide	Sux-buzone
Metho-trimepra-zine	Tamoxifen
Methox-amine	Tannic Acid
Methoxy-amine	Ten-oxicam
Methoxy-phenamine	Ter-fenadine
Methyl Salicylate	Theo-bromine
Methyl Viologene	Theo-phylline
Methylene Blue	Thiamine
Methyl-phenidate	Thi-merosal

Metricrane	Thymol
Metro-nidazole	Timolol
Tolaz-amide	Tropine
Tolbut-amide	Trypt-amine
Tolmetin	Tyrosine
Toluene	Vincamine
trans-2-Phenyl-cyclo-propyl-amine	Warfarin
Trazodone	Xylometa-zoline
Triam- cinolone	Yo-himbine
Trichlor-methiazide	Zear-alenone
Trichloro-acetic Acid Sigma Ultra	

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454X BIO-RAD Liquichek Qualitative Urine Toxicology **Positive Control** Lot No. Exp. Date :

ชื่อชุดทดสอบ **A1** = GPO Meth Cassette (1000 ng/mL) **A2** = i+Lab Methamphetamine Test Strip(1000 ng/mL) **A3** = InTec DIACNOSTIC KIT FOR METHAMPHETAMINE(500 ng/mL)

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